



PF Anti Bac Antibacterial Hand Cleaning Wipes

- Antibacterial waterless hand cleaner
- Effective against E-Coli, Leptospira (Weil's disease), etc*
- Ideal for all Utilities: gas, water, power, telecom, trains, transport, municipal, wastewater, refuse
- Heavy duty abrasive action removes even tar & bitumen
- Kind to hands – contains skin rehydrant
- Quick drying – no sticky residue
- Health & Safety – peace of mind
- Available in handy 24 wipe flat pack or 100 wipe bucket



*Tested organisms: *pseudomonas aeruginosa*, *staphylococcus aureus*, *escherichia coli*, *aspergillus niger*, *candida albicans*, *saccharomyces rouxii*.

Product Code	Description
PFI476B-ENG	PF Anti Bac Bucket 100 wipes
PFI476BBX-ENG	PF Anti Bac Box 4 x 100 wipe buckets
PFI476F-ENG	PF Anti Bac Flat Pack 24 wipes
PFI476FBX-ENG	PF Anti Bac Box 24 x 24 wipe Flat Packs

PF Anti Bac Challenge test

CFU per ml/g detected at:

Organism	0 min	30 min	60 min	3 hours	6 hours	24 hours	48 hours	7 days	14 days	21 days	28 days
<i>Pseudomonas Aeruginosa</i> NCIMB 8626	8.96 x 10 ⁶						<10	<10	<10	<10	<10
<i>Staphylococcus Aureus</i> NCTC 10788	9.12 x 10 ⁶						<10	<10	<10	<10	<10
<i>Escherichia Coli</i> NCIMB 8545	9.15 x 10 ⁶						<10	<10	<10	<10	<10
In-house contaminants	9.48 x 10 ⁶						<10	<10	<10	<10	<10
<i>Aspergillus Niger</i> IMI 149007	0.57 x 10 ⁶						<10	<10	<10	<10	<10
<i>Candida Albicans</i> NCPF 3179	8.94 x 10 ⁶						<10	<10	<10	<10	<10

PF Anti Bac independently tested for Weil's Disease prevention by Panspermia Microbiology Report CR10/1. *Leptospira interrogans icterohaemorrhagiae* Ictero1 all killed within 1 & 5 minutes @20C, in presence of 0.3 g/l Bovine albumin.